



**ROTHBAND**

PATIENT POSITIONING

## **INSTRUCTIONS FOR USE**

FOAM POSITIONING PADS – WIPEABLE COVER

VERSION 2.0

**UK  
CA CE**

## Table of Contents

<b>INTRODUCTION .....</b>	<b>3</b>
<b>INTENDED USE.....</b>	<b>3</b>
<b>REGULATORY COMPLIANCE .....</b>	<b>3</b>
<b>CAUTIONS AND WARNINGS.....</b>	<b>4</b>
<b>CLEANING INSTRUCTIONS.....</b>	<b>5</b>
1.1    DARTEX.....	5
1.2    GRIPLOCK.....	5
1.3    PVC .....	5
1.4    HEAVY DUTY NYLON.....	5
<b>DISPOSAL.....</b>	<b>5</b>
<b>CONTACT DETAILS.....</b>	<b>6</b>

## INTRODUCTION

The use, suitability and positioning of the product is beyond our control, so we do not offer guidance on these matters. The user should satisfy themselves as to suitability for intended use.

We recommend that Wipeable covered foam is always used in a protective, fluid resistant cover or, if not possible, then should be disposed of after use.

These Instructions cover the following models:

PRODUCT CODES
FPXX

## INTENDED USE

The Rothband, foam positioning pads have been designed for use with medical imaging equipment.

The positioning pads do not make any claims on maximum weight limit as it will conform to the maximum limit of the table on which it is placed.

## REGULATORY COMPLIANCE

The Wipeable foam positioning pad range complies with MDR (EU) 2017/745 (EU) and MDD 93/42/EEC (UKCA) concerning medical devices (class I) and has been designed and manufactured to comply with the following standards:

Standard	Description
EN ISO 10993-1:2009	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process (ISO 10993-1:2009)
EN ISO 10993-3:2014	Biological evaluation of medical devices - Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity (ISO 10993-3:2014)
EN ISO 14971:2012	Medical devices - Application of risk management to medical devices (ISO 14971:2007, Corrected version 2007-10-01)
EN ISO 15223-1:2016	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements.

## CAUTIONS AND WARNINGS



- The positioning aid must be x-rayed prior to use, to ensure there are no artefacts.
- Wipeable cover may cause artefact on an x-ray image.
- There are no special skills required to use this positioning aid.
- The medical professional is responsible for applying his/her best medical judgment when positioning a patient using the positioning aids.
- The positioning aids is suitable for continuous use.
- Do not place device on or near a heat source.
- The materials used in the manufacture of all components of the system comply with the required fire safety regulations.
- Intended for professional healthcare facility environments where operators with medical training are continually available when patients are present.
- Do not place heavy objects on the surface of foam pads when not in use.
- Be aware when using contrast media, as staining may cause radiographic artefact.
- The above warnings, cautions and any safety considerations should be observed on a routine and regular basis, not only upon installation.
- Patient with a high risk of pressure sore formation should not lie positioning pad for more than 10 minutes. Pressure points such as Elbows, Heels, Neck and Sacrum should be supported appropriately.

## CLEANING INSTRUCTIONS

Whilst cleaning the pads, it is always good practice to look for signs that the mattress may have been damaged. Such indicators include (but are not limited to) rips or tears and signs of fluid ingress into the interior of the pad.

### 1.1 DARTEX

Wipe with mild detergent. Wash or dry-clean using detergent up to 95oC. (remove foam element first) For disinfection use sodium hypochlorite up to concentrations of 10,000 parts per million of active chlorine. Do not use biological or phenolic-based cleaners. Wash and dry completely before storage. Some surface wrinkling may occur.

### 1.2 GRIPLOCK

Wipe with mild detergent. Wash or dry-clean using detergent up to 75oC. (remove foam element first) For disinfection use sodium hypochlorite at no more than 1000 parts per million of active chlorine. Do not use biological or phenolic-based cleaners. Wash and dry completely before storage. Some surface wrinkling may occur.

### 1.3 PVC

Wipe with mild detergent. Avoid chlorinated cleaners and it attacks the plasticisers and can result in a hardened coating. For disinfection use Rothband apron cleaner. Do not use biological or phenolic-based cleaners. Wash and dry completely before reuse. Some surface wrinkling may occur.

### 1.4 HEAVY DUTY NYLON

Wipe with mild detergent. For disinfection use sodium hypochlorite at no more than 1000 parts per million of active chlorine. Do not use biological or phenolic-based cleaners. Wash and dry completely before storage. Some surface wrinkling may occur.

## DISPOSAL

We recommend final disposal to be as clinical waste.

## CONTACT DETAILS



WSR Medical Solutions  
(Trading as Rothband)  
Unit 2 Hargreaves Mill  
Hargreaves Street  
Haslingden  
Lancashire  
BB4 5RQ  
UK



Eurolink Europe Compliance Limited, Dublin  
D02 AY86

Tel: 01706 830086  
e-mail: [sales@rothband.com](mailto:sales@rothband.com)