

# **MOBILE AND THEATRES**

## MOBILE AND THEATRES 37





## **SLIDERERS - DISPOSABLE CASSETTE COVERS**

Made from an extremely strong but slippy plastic which both Assists in cassette placement & in Infection Control. Easier placement & removal of cassette/plate. Reduces back, arm & neck strain for staff. Minimises movement required by patient. Protection of cassette from fluids & contamination. Protection of staff AND patients. Dispose of as clinical waste. Rail available for wall / Mobile machine mounting

ORDER CODE	Description
SLIDA	Sliderers 0-100
SLIDB	Sliderers 150-1000
SLIDC	Sliderers 1050+ (state volume)

# **PAG - DR PLATE HOLDER**

The Protect-A-Grid DR panel premium grid encasement allows for safe & effective examinations when using Direct Radiographic plates.

The PAG-DRP is a drop-on style that snaps on & off your existing DR panel quickly & effectively.

Please provide us with details of your equipment manufacturer & the size of the plate you require.

ORDER CODE	Description
PAG	Encasement without grid
PAG2430AA	PAG with grid, 24 x 30, 8:1/40 lines AL/AL
PAG3540AA	PAG with grid, 35 x 43, 8:1/40 lines AL/AL
PAG3540CA	PAG with grid, 35 x 43, 8:1/40 lines CF/AL
PAG3570AA	PAG with grid, 35 x 43, 8:1/70 lines AL/AL
PAG3570CA	PAG with grid, 35 x 43, 8:1/70 lines CF/AL
PAGDRX1-A	PAG designed for DRX1, 215/10:1
PAGDRX1-B	PAG designed for DRX1, 103/10:1





# PADDED CASSETTE BACKREST

For chest examinations. Use behind patients when on trolleys/beds. Helps obtain an erect position. Fits any DR Plate, please specify. Wipeable & removable Nylon outer covers. 2 handles for easy carrying/positioning.

Approx 150 angle 40 x 50 x 15cms

ORDER CODE PCBR





## **MOBILE PARKING SPACE**

Ideal for identifying the parking space of mobile xray equipment to ensure access to charging points. Hardwearing textured surface Suitable for almost all indoor floor materials R12 slip resistance & B1 fire safety certified Fully cusomtisable Easy to apply

ORDER CODE MOBPS

## MOBILE AND THEATRES 39

C-ARM MONITOR COVER REMOVE BEFORE USE FOR X-RAY USE ONLY

> C-ARM INTENSIFIER COVER REMOVE BEFORE USE FOR X-RAY USE ONLY

ARM CONTROL PANEL COVI REMOVE BEFORE USE FOR X-RAY USE ONLY

# **C-ARM COVER SET**

Designed to offer scratch protection to imaging equipment whilst in storage on busy operating theatre corridors. The cover contains a tough polyprop board with cushioning.

Includes:

- Intensifier Cover
- Monitor Cover
- Control Panel Cover

Please contact your representative for help in measuring the required dimensions

ORDER CODE CARMCSET





