

MR COMPATIBILITY CERTIFICATE

The following product has been tested and found to be MR Safe in a 3 Tesla Environment and lower.

03-WCA26DFB-650 - MRI Non-Magnetic Bariatric Wheelchair 650 lbs Capacity

03-WCB26DFB-850 - MRI Non-Magnetic Bariatric Wheelchair 850 lbs Capacity

The product is made of Non-Magnetic Stainless Steel, Brass and Plastics.

Results: Non-clinical testing demonstrated that each product is MRI Conditional and can be used in the MRI environment according to the following conditions:

- Static magnetic field of 3-Tesla or less
- Spatial gradient magnetic field of 720-Gauss/cm or less

There are thousands of these units in use in 3 Tesla MRI Suites around the world and used on a daily basis.

These can be used in the MRI Magnet Room with the following Conditions:

This product is intended for use inside of the MRI environment (e.g., in the MR system room, close to the scanner). However, it should not be utilized directly inside of the MR system (e.g., inside of the bore of the scanner), itself.

Thank you for the interest in our products.

Eric B Schatzlein

General Manager
Magmedix Inc

