

US-11

# Scrotal Ultrasound Phantom

## Master techniques for a safe and painless examination

Hands-on training to accurately diagnose scrotum abnormalities

### Product Supervision


Charlotte Henningsen, MS, RT(R), RDMS, RVT,  
FSDMS, FAIUM  
Chair & Professor - Sonography Department  
Adventist University of Health Sciences

### FEATURES

- 1 | Two types of interchangeable scrotums: normal and pathology
- 2 | Excellent ultrasound image quality

### SKILLS




- | Patient positioning and preparation for examination
- | Scrotal ultrasound screening
- | Visualization of testicular cancer



## KEY FEATURES




### Normal phantom

- scrotum, testicle, epididymis and penis



### Pathological phantom

-scrotum, testicle, epididymis and penis  
-testicular cancer (each one in left and right testis, 10mm dia.)



## DESCRIPTIONS

### SPECIFICATIONS

Size: 34x33x24cm / 13.4x13x9.4 inch Weight: 3.6kg / 7.9lbs

### MATERIALS

- 1 pelvic body
- 1 normal scrotum
- 1 pathological scrotum
- 1 storage case
- 1 instruction manual

### REPLACEMENT PARTS

- 41911-010 1 normal scrotum unit
- 41911-020 1 pathological scrotum unit

\* Ultrasound scanner is not included.

